

저주파 자극기

TRANSCUTANEOUS ELECTRIC NERVE STIMULATION (TENS) —

Biomix HSP-601

HSP-601은 저주파 자극을 이용하여 통증완화에 사용하는 기기로 흡입도자를 채택하여 부드럽고 6개 채널에서 동시에 출력이 가능하여 효과적인 치료가 가능합니다.

Through the operation of switch, it is possible to control the wave type of electric current with low frequency, the operating time, the pressure strength of absorbing motor, the start/the end of absorbing motor, and the start/the end of output. Also, it is possible to carry out from six channels at the same time.


Feature

- Fasteners are unnecessary due to suction electrode.
- Use 6 channels at the same time.
- Various therapy area with small and large electrode.
- Low frequency is used for stimulation to relieve pain.

Specifications

Power	AC 220V, 50/60Hz
Power Consumption	40W
Safety, Isolation	Class I, Type BF
Dimension	430x315x760(H)mm
Weight	20Kg

Accessory



Suction cup

Various therapy area and easy-fixing with small and large suction electrode.



HSP-601