

간섭 전류형 저주파 자극기

INTERFERENTIAL CURRENT THERAPY (ICT)

BM-750 / BM-760

근육통 완화 등의 목적으로 피부표면에 부착한 전극을 통하여 피부저항이 거의 없는 서로 다른 중주파를 가하는 간섭전류형 저주파 자극기로 발진부, 출력조절부, 증폭부, 타이머, 지시부, 전원부, 안전장치 및 자극 도자로 구성되어 있으며 좌우측에 별개의 출력단자가 있어 2인 동시 자극이 가능합니다. 5가지의 자동 프로그램과 25가지의 수동 모드로 편리하고 선택적 치료가 가능합니다.

This device is for "Interferential Current Therapy" that helps to relieve and reduce pain, such as muscular pain, using various medium frequencies coming from electrodes attached on skin surface. BM-760 is composed of generator, output controller, timer, indicator, power supply, safety device, and stimulus electrodes. The device may use for two patients at the same time with two isolated outputs located in left and right sides and has 5 built-in auto programs and 25 manual modes.

Feature

- 5가지의 자동 프로그램으로 사용이 편리
Easy operation with 5 auto programs.
- 25가지의 수동 모드로 환자의 선택적 치료 가능
Effective therapy for each patient with 25 manual modes.
- 2인 동시 치료 가능
Treatment for 2 patients at the same time.
- 2채널이 분리되어 있어 각기 다른 부위에 사용 가능
Two isolated outputs for individual treatment.

Specifications

Power	AC 220V, 50/60Hz
Power Consumption	200VA
Safety, Isolation	Class I, Type BF
Output Wave	Sine Wave
Max Output Current	80mA
Stimulus Frequency	0~1100Hz
Basic Wave Frequency	3KHz, 4KHz, 5KHz
Dimension	338x458x950(H)mm
Weight	26Kg

* BM-760 (L: 6-pole R: 6-pole), BM-750 (L:4-pole R: 6-pole)

Accessory



흡인형 도자
(Electrode Suction Cups with cables)



저주파 패드 별도구매
(Electrode Suction Cups with cables
Optional)



BM-760