

저주파 자극기

TRANSCUTANEOUS ELECTRIC NERVE STIMULATION (TENS)

H-3000

근육통 완화 등을 목적으로 피부 표면에 부착한 전극을 이용하여 저주파전류를 가하는 저주파 치료기로서, 2개의 채널을 2인이 동시에 사용하도록 별도의 출력 조절이 가능한 보급형 기기입니다. 다양한 치료 프로그램을 선택하여 치료 효과를 극대화시켰습니다.

H-3000 is a TENS machine to use for two patients individually. Electrode suction cups and Electrode pads can be applied separately at the same time.

Feature

- 패드와 흡인컵 각각의 독립적 사용이 가능한 2인용 치료기(흡인형)  
The pad and suction electrode cup could be individually used at the same time.
- 6개의 자동프로그램으로 부위별 치료가 가능  
Treatment for main part is possible with 6 auto programs.
- 3개의 고정프로그램으로 치료출력을 가변하여 부드러운 치료가 가능  
Soft treatment with 3 fixed frequency programs by changing the output power.
- 25종류의 주요한 치료부위 프로그램과 증상에 따른 프로그램 선택으로 질환별 치료효과 극대화  
Treatment for special part with 5 part programs and for symptoms with 5 symptom programs.
- 사용자 편의를 위한 모든 기능을 쉽게 확인할 수 있는 표시부  
Optical Display for all functions and rack to store suction cups/electrodes for users' convenience.


Specifications

Power	AC 220V, 60Hz
Power Consumption	40W
Safety, Isolation	Class I, Type BF
Output Current	40mA
Vacuum Pressure	max 250mmHg
Timer	0~15min
Dimension	450x350x160(H)mm, 13Kg

Accessory

흡인도자 2세트  
(4 Electrode suction cups with cable, sponge)

패드도자 2세트  
(4 Electrode pads with cable, band)



저주파패드  
(Electrode Pads with cables)



흡인형도자  
(Electrode Suction Cups with cables)



H-3000  
CFDA