

전동식 정형용 견인장치

ORTHOPEDICS TRACTION UNIT

HT-101, HT-102 Biotrack

물리치료에서 빼놓을 수 없는 견인요법.

무엇보다도 중요한 것은 안전성이며 매일 사용해야 하므로 부드럽고 사용하기 쉬워야 합니다.

HT-101(1인용)과 HT-102(2인용)은 안전하고 쉽게 사용할 수 있도록 설계되었고, 전동식 장치를 이용하여 경추견인과 요추견인을 동시에 할 수 있으며 견인의 상승과 하강속도, 환자의 체중, 작동시간을 조절할 수 있도록 업그레이드 되었습니다.

온열 기능이 있고 안전스위치가 제공됩니다.

The electrically-powered orthopedics traction unit is designed to use softly and easily.

This appliance can track the cervical and the waist vertebrae individually at the same time. HT-101 is for one patient and HT-102 for two patients. An Emergency stop switch for a patient is provided.

Feature

- 4 traction modes · Computerized control.
- Hot-pack effect in the back area of the treatment table.
- Convenient operation with shift lock.
- If a patient moves during treatment, the device adjusts automatically to the Patient's position.
- It is always maintaining fixed pulling power for treatment.


Specifications

Power	AC 220V, 60Hz
Power Consumption	120W(HT-101)/200W(HT-102)
Safety, Isolation	Class I, Type B
Pulling Power	1~80Kg
Pulling Speed	7mm/sec
Timer	1~99min
Hold/Standby Time	1~99sec
Thermal effect	{Level 1} 40°C / {Level 2} 60°C
Dimension	2425 x 580 x 2640(H)mm
Weight	60kg(Control Unit), 65kg(Bed)


* HT-102S (with one bed), HT-102D (with two beds)



HT-101
Single system



HT-102
Dual system



HT-101



HT-101 Biotrack
Single Control Panel



HT-102 Biotrack
Dual Control Panel

HT-102
CFDA

Accessory

